



January 2022

Homechoice Claria with Updated Software Shipping Soon

Dear Valued Patient,

In the coming months the Homechoice Claria cycler software will be updated to version 12.002 in Baxter shipment centers. If you receive a replacement cycler you may notice that the way alarms are displayed on the screen has changed. Some alarms are now called “Notices” and new symbols will appear on the screen when an alarm is active and muted. This update does not require any action on your part.

Thank you for trusting Baxter to support you on your PD journey.

Included Product:

Product Code	Description	Serial Numbers
5C6M40	Homechoice Claria APD system	All

Alarm Active Symbol: Flashing Triangle



Alarm Muted Symbol: Crossed-out Bell



Sincerely,

Alexandra Duffy
PD Products
Baxter Healthcare Corporation

Homechoice Claria APD system is intended for automated dialysis solution exchanges in patients with renal failure undergoing peritoneal dialysis, at home. **Rx Only.** For the safe and proper use of the devices mentioned herein, refer to the Instructions for Use.

Baxter and Homechoice Claria are trademarks of Baxter International Inc.